



Recent Results for the µPIC Neutron Imaging Detector

KYOTO UNIVERSITY, COSMIC RAY GROUP

Joe Parker, K. Hattori, S. Iwaki, S. Kabuki, Y. Kishimoto, H. Kubo, S. Kurosawa, Y. Matsuoka, K. Miuchi, H. Nishimura, T. Sawano, A. Takada, T. Tanimori, K. Ueno

JAPAN ATOMIC ENERGY AGENCY, MATERIALS AND LIFE SCIENCE FACILITY DIVISION M. Harada, T. Oku, T. Shinohara, J. Suzuki

> HOKKAIDO UNIVERSITY T. Kamiyama, K. Kino, Y. Kiyanagi

Neutron imaging detector prototype (µNID)



- * TPC measures 3D proton-triton tracks.
- Time-over-threshold to estimate energy deposition.
- Gas gain < 1000 for neutron imaging.</p>
- Efficiency up to ~30%, spatial res. of ~120 μm, time res. of ~1 μs.
- Detector remains operable for over 1 year on single gas filling.

FPGA Encoder with time-overthreshold



- Fast 3D tracking and 'energy' measurement
- Improved position accuracy and background rejection.

Spatial resolution with template fit

- Proton range ~3 times triton.
- Proton-triton separation, position correction essential for spatial resolution.
- Proton direction, neutron position determined by fit.

Measured TOT distribution



Templates generated with GEANT4 simulation of detector.

Data taken at NOBORU in Feb. 2011.



- Improvement possible (10~15%) by optimizing gas for reduced diffusion/shorter track-lengths.
- Recent data (NOBORU, March 2012) suggests some improvement (~5%) with lower gas gain.

Neutron-gamma separation

- Gamma rejection studied using 1-MBq ¹³⁷Cs source.
- Data taken over 24 hours at a gas gain of ~600.

Pulse-width sum after track-length cut



Fraction of detected γ 's surviving neutron cuts < 10⁻⁶ (effective gamma sensitivity of < 10⁻⁹ at gain of ~600).

Sensitivity of $< 10^{-12}$ can be achieved by reducing gain by half without loss of neutron efficiency.

- * Both neutrons and gammas are detected (γ efficiency < 10⁻³).
- Separation by total time-over-threshold and 3D track length.

n and γ event rates (¹³⁷Cs + ²⁵²Cf)



Recent measurements

* Resonance absorption.
* Bragg-edge transmission.
* CT test measurement.



Neutron resonance absorption



- * Exposure: 5 hrs.
- **★** TOF: 0 2.2 ms.
- ***** ~10 kcps.
- ***** Live time: 60%.
- Measurement of neutron TOF allows selective imaging of nuclides via resonance absorption.
- Good time resolution is essential.

Image vs. TOF (0 ~ 1 ms)





Data taken at NOBORU in March 2012.

Bragg-edge transmission

2D image of katana



Region used for Bragg-edge study

- **FOF** (ms) **TOF projection for** one slice 0.75 0.7 0.65 0.6 0.55 0.5 3.5 4 4.5 5 5.5 6 6.5 7 7.5 TOF (ms) Edge fit: error function with exponentials above and below. Simplified version of function from J.R. Santisteban et al., J Appl Cryst 34 (2001).
- Fragment of a katana (15th or 16th century).
- Bragg edge positions indicate some difference between cutting edge and spine.

Data taken at Hokkaido University in February 2012.



CT test measurement





Improvements

* New ASICs and encoder for faster, more compact DAQ.

DAQ upgrade

Current DAQ setup



- * ASD racks, cables, and encoder module.
- Output through single FPGA limits DAQ rate to ~10 MHz (neutron rate of ~150 kHz).
- Slow VME-PC transfer limits live time (15 ~ 60%).

New encoder module



118mm × 220mm

128 channels per board

- * ASICs and FPGA on single board.
- * Four boards replace ASDs, encoder, cables.
- DAQ rate increase by factor of 5 or more.
- * Will begin testing with neutron system soon.

at Kyoto Univ.

cm test system

VME-PC transfer

- Implement faster VME readout (block transfer, double-buffer).
- Block transfer gives 1.5X increase; double-buffer not yet tested.
- * For future, switch to SiTCP.

Summary

µPIC-based time-resolved neutron imaging detector.

 For radiography, neutron resonance absorption imaging/spectroscopy, Bragg-edge transmission, CT imaging, SANS.

Detector performance.

- * Spatial resolution of $105 \sim 130 \mu m$; time resolution of $\sim 1 \mu s$.
- Very small effective gamma sensitivity of < 10⁻¹².
- After upgrading DAQ and optimizing gas, maximum neutron rate up to 1 Mcps and spatial resolution < 100 μm.

* Next

- Develop and test new FPGA code for neutron imaging.
- Possible beam test of new DAQ at NOBORU this winter.
- Consider design of next version of the detector.

Extra slides



Resonance absorption: Ag-In-Cd alloy



- Plate thickness: 3 mm.
- * Exposure: 2 hrs.
- * TOF gate: 0 3 ms.
- * Neutron rate: ~10 kcps.
- * DAQ live time: 70%.

Image of ASTM indicator

Image taken with µPIC (100 µm bins)



- * Exposure: 3 hrs.
- * No TOF gate.
- * ~120 kcps.
- * Live time: 14%.



Data taken at NOBORU in March 2012.

X-ray provided with sample





DAQ and FPGA logic

Amplifier-Shaper-Discriminators (ATLAS, KEK)



DATA ENCODING

- Two words per pulse.
- * 'edge bit' saved with each data word.



PROTON-TRITON TRACKS



- Simultaneous measurement of position and 'energy deposit' at high rates.
- Excellent background rejection capability.

Neutron efficiency

- Neutron efficiency as a function of neutron energy.
- Determined from GEANT4 simulation.
- Loss of peak efficiency due to large dead layer in current configuration (2.5-cm).



Drift cage configuration	Drift height	Dead layer	Efficiency at 25.3 meV	Peak efficiency
5-cm	5.0 cm	0.8 cm	0.35	0.75 (0.7meV)
2.5-cm	2.5 cm	3.3 cm	0.18	0.27 (3meV)
2.5-cm(2)	2.5 cm	0.8 cm	0.20	0.65 (0.35meV)

Optimize gas for improved position resolution

Strategies for improving position resolution

- Shorten proton-triton track lengths.
- Reduce diffusion of drift electrons.

Gas parameters determined by MAGBOLTZ. Resolutions estimated with GEANT4.

	Pressure (atm)	Drift velocity (µm/ns)	Transverse diffusion (µm/cm ^{1/2})	p-t track length (mm)	Expected improvement in resolution
Ar:C ₂ H ₆ : ³ He (63:7:30)	2	23.1	273	7.9	(114 µm)
Ar:C ₂ H ₆ : ³ He (63:7:30)	3	23.4	231	5.3	~15%
Xe:C ₂ H ₆ : ³ He (50:20:30)	2	29.4	183	5.0	~15%
Ar:CO ₂ : ³ He (50:20:30)	2	22.5	107	7.4	~15%

Preliminary tests with Ar:CO₂ mixtures show relatively low gain.

Also considering other CO₂-based gas mixtures.

Gas study – Ar:CO₂-based mixtures

Compared following mixtures at 1 and 2 atm:

Ar:C₂H₆:He (63:7:30)

Ar:CO₂:He (50:15:35)

Ar:CO₂:C₂H₆:He (50:10:5:35)

- C₂H₆ included to increase gain and stability.
- ^{*} ⁴He used in place of ³He.

CO2-based mixtures require higher anode voltages.

Next steps:

- Check gain variation.
- Add small amount of ³He and measure proton-triton tracks.



Gain variation and pixel pitch

- * Normal gain variation of μ PIC: 4%(σ).
- * Gain variation of prototype neutron detector: $16\%(\sigma)$.
- Gain variation contributes 10~20% to resolution (GEANT4 simulation).
- Recent improvements in manufacturing procedure may reduce intrinsic gain variation of µPIC.



- Change arrangement of pixels (physical dimensions unchanged).
- Moderate reductions in pixel pitch produce corresponding reduction in position resolution (GEANT4).



Previous experiments

SMALL-ANGLE NEUTRON SCATTERING (2009) Spherical SiO₂ nanoparticles (ϕ 200µm).



RESONANCE IMAGING (2009) Square cobalt sample.



BRAGG-EDGE TRANSMISSION (2010) Welded steel plate.



Resonance absorption

- Transmission for sheets of Mn, Mo, Ag, Ta, and In (thicknesses 10µm ~ 1mm).
- Time resolution of ~0.6 µs and good background rejection can distinguish peaks near beginning of pulse.



Bragg-edge transmission



- Transmission for Fe powder (>99% pure, grain size < 325 μm).
- Sample thickness of 1.6 cm.
- Bragg-edges are clearly visible.
- Edge spacing is consistent with expected BCC crystal structure.

Data taken at NOBORU, J-PARC in Feb. 2011.