

Neutron Imaging Detector Based on the μ PIC Micro-Pixel Chamber


Kyoto University, Cosmic Ray Group

Joe Parker, K. Hattori, S. Iwaki, S. Kabuki, Y. Kishimoto, H. Kubo, S. Kurosawa, Y. Matsuoka, K. Miuchi, H. Nishimura, T. Sawano, A. Takada, T. Tanimori, K. Ueno

Japan Atomic Energy Agency, Materials and Life Science Facility Division



M. Harada, T. Oku, T. Shinohara, J. Suzuki

Time-resolved neutron imaging detector




- * Ar:C₂H₆:³He gas mixture at 2 atm.
- * Gas gain < 1000 for neutron imaging.
- * Efficiency of 18% at 25.3 meV and up to 65% at 0.35 meV (2.5-cm gas depth).
- * ~105 μm spatial resolution; 0.6 μs time resolution; <10⁻¹² gamma sensitivity.
- * 2 years operation or more on single gas filling.

FPGA encoder with time-over-threshold (TOT)




Improved position accuracy and background rejection via detailed tracking.

DAQ and FPGA logic




Rate test at NOBORU, J-PARC (2012)

Neutron intensity varied by adjusting size of slit opening



Data encoding

- * Two words per pulse.
- * 'Edge bit' saved with each data word.




DAQ shows good linearity up to a data rate of 10~12MHz (neutron rate of ~200kcps).

Deviation due to internal FIFO overflows in FPGAs.

Fast neutron and gamma rejection

TOT distribution for scattered proton

Fast neutrons from ^{252}Cf source




Track-length and energy dep. can be similar to neutron events.

Distinguished by shape of TOT distribution.

Event-by-event TOT sum (1-MBq ^{137}Cs)



Gas gain ~ 600 , track-length cut applied

Gamma efficiency $< 10^{-3}$ at gain of 600.



Fraction of detected γ 's surviving neutron cuts $< 10^{-6}$ (effective gamma sensitivity of $< 10^{-9}$ at gain of ~ 600).

TOT sum with no cuts and n and γ event rates ($^{137}\text{Cs} + ^{252}\text{Cf}$)




Sensitivity of order 10^{-12} or less can be achieved at reduced gain without loss of neutron efficiency.

Spatial resolution


- * Proton-triton separation essential for spatial resolution.
- * Proton direction, neutron position determined by 'template' fit.

Image taken at NOBORU, J-PARC in Feb. 2011.
Gas gain ~470, Exposure time 29 mins.




Position from mid-point of track.
Resolution: ~1mm(σ)

Measured TOT distribution



Bin size: 200 μm \times 200 μm .

+PID



Bin size: 80 μm \times 80 μm .

Preliminary!
Resolution with PID:
104 to 124 μm (σ)
(Final resolution depends on tracking cuts.)

$$TOT(x, m) = A \text{ Template}_m [\hat{x} f_w(x - x_N)]$$


- * Templates generated with GEANT4 simulation of detector.

Spatial resolution

Images taken at NOBORU, J-PARC in March 2012.

Resolution from knife-edge test


Cadmium test chart, gas gain ~470



	Res(μm)	Eff(%)
A	131 ± 3	100
B	125 ± 3	66
C	111 ± 4	32



- A. Gamma-rejection cuts.
- B. + Z-consistency.
- C. + Two peaks in TOT.

Resolution vs gas gain





Improves by $\sim 3\mu\text{m}$ per 100 decrease in gain.

Modulation transfer function



Uniformity of resolution

Evaluated over 3cm x 3cm area



After correcting for beam divergence

Mean: $104\mu\text{m}$
RMS: 5.3%

Improvement of spatial resolution


Gas properties

- * Shorten proton-triton track lengths.
- * Reduce electron diffusion.

Improvements to μ PIC

- * Improve gain variation (better manufacturing controls).
- * Reduce pixel pitch.

Measured gain variation of prototype



Simulation results for different gas mixtures

	Pressure (atm)	Transverse diffusion ($\mu\text{m}/\text{cm}^{1/2}$)	p-t track length (mm)	Expected improvement in resolution
Ar:C ₂ H ₆ : ³ He (63:7:30)	2	273	7.9	(105 μm ~)
Ar:C ₂ H ₆ : ³ He (63:7:30)	3	231	5.3	~15%
Xe:C ₂ H ₆ : ³ He (50:20:30)	2	183	5.0	~15%
Ar:CO ₂ : ³ He (50:20:30)	2	107	7.4	~15%

Gas parameters determined by MAGBOLTZ.
Resolutions estimated with GEANT4.

- * Reducing gain variation improves spatial resolution by 5~10% and uniformity by 15~30%.
- * In low gain operation, pixel pitch can be reduced to ~200 μm .
- * GEANT4 studies are ongoing.


Resolution after diffusion and gain variation improvements: 80~90 μm .

Plus increased pixel pitch: < 70 μm .

Preliminary!

DAQ upgrade

Current DAQ



Many cables transfer signals to FPGA encoder


Output through single FPGA limits DAQ rate to 10~12MHz (~200kcps)

VME-PC transfer limits live-time

VME-PC transfer


- * Implement faster VME readout (block transfer, double-buffer).
- * Block transfer gives 1.5X increase; double-buffer not yet tested.
- * For future, switch to SiTCP.

New encoder module




- * 128 channels per board.
- * ASICs and FPGA on single board.

Expect rate increase by factor of 4 or more.



Cosmic-ray muon event




TPC mode (digital hit data + FADC)

Recent measurements

- * Resonance absorption.
- * CT test measurement.



Neutron resonance absorption




Exposure: 5 hrs
Rate: ~10 kcps
Live time: 60%


External TOF
gate applied (0 –
2.2ms)

- * Measurement of neutron TOF allows selective imaging of nuclides via resonance absorption.
- * Good time resolution is essential.


Image vs. TOF (0 ~ 1 ms)





Transmission in sample region




CT test measurement



Goniometer with θ -rotation from -155° to 155° .




- * 9 angles (6 independent), 54 minutes/angle.
- * TOF: 0 – 3 ms
- * Neutron rate ~ 10 kcps, Live time 60%.
- * CT reconstruction using simple back-projection method.




Summary

- * μ PIC-based time-resolved neutron imaging detector.
 - * For radiography, neutron resonance absorption imaging/spectroscopy, Bragg-edge transmission, CT imaging, SANS.
- * Detector performance.
 - * Spatial resolution of 105 μm (100% reconstruction eff.) to 120 μm (32% reconstruction eff.).
 - * Time resolution $\sim 0.6 \mu\text{s}$; very small effective gamma sensitivity $< 10^{-12}$.
 - * Can operate for ~ 2 years on a single gas filling.
 - * After upgrading DAQ and optimizing gas, maximum neutron rate on order of 1 Mcps and spatial resolution $< 100 \mu\text{m}$.
- * Next
 - * Develop and test new FPGA code for neutron imaging.
 - * Beam test of new DAQ at NOBORU this winter.

Extra slides



Resonance absorption: Ag-In-Cd alloy



- * Plate thickness: 3 mm.
- * Exposure: 2 hrs.
- * TOF gate: 0 – 3 ms.
- * Neutron rate: ~10 kcps.
- * DAQ live time: 70%.





Image of ASTM indicator

Image taken with μ PIC (100 μ m bins)



- * Exposure: 3 hrs.
- * No TOF gate.
- * ~120 kcps.
- * Live time: 14%.

X-ray provided with sample

